

## Export a design

---

### 1 Why export?

- This is used to export a study design to a CDISC ODM file, with or without CDISC SDM and Viedoc extensions. It can be used for import in another project or another instance of Viedoc, e.g. a training instance. As the file is CDISC compliant, it can also be used in other systems equally compliant with CDISC standards.

A screenshot of the 'Design Settings' dialog box. It has two tabs: 'Details' and 'Export Design'. Under 'Export Design', there are three options: 'Include Viedoc extensions?' with a checked checkbox, 'Include SDM extensions?' with a checked checkbox, and 'Export' with a red arrow pointing to it.

Design Settings

Details Export Design

Include Viedoc extensions?

Include SDM extensions?

Export

CDISC SDM contains study workflow information and Viedoc extensions are Viedoc specific settings that cannot be described as part of the CDISC standards. If Viedoc is the target system, both checkboxes should always be checked.