


Export a design

1 Why export?

- This is used to export a study design to a CDISC ODM file, with or without CDISC SDM and Viedoc extensions. It can be used for import in another project or another instance of Viedoc, e.g. a training instance. As the file is CDISC compliant, it can also be used in other systems equally compliant with CDISC standards.



The screenshot shows the 'Design Settings' dialog box with two tabs: 'Details' and 'Export Design'. The 'Export Design' tab is active. It contains two checkboxes, both of which are checked: 'Include Viedoc extensions?' and 'Include SDM extensions?'. Below these checkboxes is an 'Export' button, which is highlighted by a red arrow.

CDISC SDM contains study workflow information and Viedoc extensions are Viedoc specific settings that cannot be described as part of the CDISC standards. If Viedoc is the target system, both checkboxes should always be checked.