

# 互 Data mapping for import of data

# About importing data into Viedoc 4

- Viedoc offers the possibility to import data, for example laboratory data, via the Viedoc 4 Data Import Application. The Viedoc 4 Data Import application does the following:
  - it converts the provided data into CDISC ODM clinical data format using a data mapping file, and
  - it pushes the data into Viedoc through the Viedoc 4 API.
- In order to import data into Viedoc, the following is needed:
  - a data mapping file
  - a configuration file
  - the Viedoc Data Import application
  - a valid Viedoc username and password, and a study-specific Viedoc 4 API client key
  - the data to be imported, supplied in a CSV (comma-separated value)-file
- The data mapping file defines how the external data will be mapped to form items in Viedoc. Internally, the data mapping
  will be stored in CDISC Define-XML format.
- The configuration file is an xml file that defines:
  - where to find the data mapping file,
  - where to find the data that will be imported,
  - into which Viedoc study data will be imported,
  - which API instance the data should be imported into,
  - the login credentials that will be used when importing the data.

In the configuration file, you can also define whether you would like new subjects to be created automatically during the data import, when the imported data contain data for a subject that has not been added to the study yet. It is also possible to define whether you would like visits to be initiated during the data import, when the imported data contain data of visits that have not been initiated yet.

One configuration file can contain the import configurations for multiple studies.

- The import of laboratory data into Viedoc using the Viedoc 4 Data Import Application involves the following steps:
  - 1. Creating a data mapping file in Viedoc Designer.
  - 2. Creating a Viedoc 4 API client key in Viedoc Admin.
  - 3. Creating a configuration file (not in Viedoc).
  - 4. Preparing the work folder.
  - 5. Installing the Viedoc 4 Data Import Application.
  - 6. Dropping data into the work folder.
  - 7. Running the Viedoc 4 Data Import Application.

#### More information

- This lesson describes how to create a data mapping file in Viedoc Designer.
- For more information on how to use the Viedoc 4 Data Import Application, see the Viedoc 4 Data Import Application document. The Viedoc 4 Data Import Application document is available on request, please contact PCG Solutions AB.
- Instructions on how to create a data mapping file and how to import data using the Viedoc 4 Data Import Application can also be found in our **video tutorial**. Click <u>here</u> to watch the video.

## 3 Description of the data mapping file

- The data mapping file is created in Global design settings in Viedoc Designer. It describes each column in the CSV file, and how and where the data should be imported into Viedoc.
- A separate data mapping file should be created for each study. It is possible to create more than one data mapping per study.
- The data mapping window in Viedoc Designer has the following main fields (see image):
  - 1. Data Mapping Name, a name (free text) for the data mapping.
  - 2. **Domain Name**, the SAS domain that will be stored in the Define-XML file according to the CDISC standard.
  - 3. Data mapping table, that has two main parts:
    Imported file structure, describing each column of the CSV file to be imported.
    - **Viedoc**, describing how and where the data should be imported into Viedoc.
- The data mapping table has the following columns:

Data Happing Name Data Happing Name						Donain name Safety Laboratory Parameters							
lespe	orted file structure						Vieda	4					
	Column name		Description		Link to		14	Destination (Viedoc expression)			cı.		
1	90		Subject key										
-1			Subject key					Eubjectioyi					
-2			Country site and sequence					ICountryCodel-(SiteCode)-(SiteSubjectSeqNo)					
8	HaemCollected	1			Nove			ISTHIS LOLDHARYN)					
	CL1: yes							04.1		8			
	CL2 ro							CL2 2		8			
3	HermDate		Date and time sample collection		Nove			(STHIS LB LBHAEDT)					
4	HaemResult.				Note			(STHIS LB LBHADRES)	1				
	CL1: normal							CU: 1		8			
	CL2 abnormalizaignif							CL2 2		8			
	CL3 abnormabignit							CL3 3		8			
5	HaemAbrionnal				None			ISTHIS USUBHARSPECT					
	CL1: hemoglabin							CUL 1		8			
	CL2 hematocrit							C15 5		8			
	CL3 MCV							C18 3		8			
	CL4 W9C							CL4 4		8			
	CLS: planelets							CLS 5	1	8			

Column name	Description
#	Number of the column in the CSV file to be imported.
Column name	Name of the column in the CSV file to be imported.
Description	Description of the parameter in that specific column to be imported (free text).
Link to	Links the content to a parameter defined in another column in the CSV file. This is mainly used when importing data in a tall-skinny format.
IM	Item mapping: inserts more rows into the table so that one CSV column can be mapped to multiple destinations in Viedoc.
Destination	Address in Viedoc into which the data should be imported. This directs the data to the correct patient, visit, form and field.
CL	Code list: a list of codes and their definition.

# 4 Creating the data mapping file

- Follow the steps below to create a data mapping file.
  - 1. In Designer, scroll to the study for which you would like to create a data mapping file.
  - 2. In the **Global design settings** field, click **Edit** to open the global design settings.
  - 3. In the Data mappings field, click Edit to open the data mappings overview.
  - Click Add new data mapping, or if a data mapping has already been created, select the data mapping that you would like to edit.
  - 5. Type a name for the data mapping in the **Data Mapping Name** field.
  - 6. Type a domain name in the **Domain name** field.
  - 7. Click **Add new row** to add a new row to the table, and fill out:
    - Column name of the first column of the CSV file.
      - Description of the first column of the CSV file.
      - Destination in Viedoc of the data in the first column of the CSV file.

For more information on how to correctly map the patient ID, visit ID and laboratory data, refer to the sections below.

8. Repeat step 7 with the next column in the CSV file, until all the columns in the CSV file are described in the data mapping table.

If you would like to link the contents of the current row to the contents of another row, select the row to link to from the **Link to** dropdown menu. This is mainly used for data in tall-skinny format. Choosing **Link to** automatically creates rows with the same code list items as the row that is being linked to. Linked rows can be updated by clicking on the refresh button.

- 9. Click **Save changes**, and click **Close** to close the data mapping table.
- 10. Click **Close** to exit the data mappings overview.

- 11. Click **Publish settings** in the global design setting window to publish the changes.
- 12. In the Data mappings field, click Edit to open the data mappings overview.
- 13. Click the **Download** icon behind the data mapping that you just created. An xml file will be downloaded that contains the data mapping.
- 14. Save the xml file in the work folder.
- You can enter more mapping details by clicking the +-icon in the **IM** column. This adds extra rows to the table that can be used to map the data to more than one destination in Viedoc.
- You can enter code list items by clicking the +-icon in the CL column. Code list items define how values in the CSV file should be translated into values in Viedoc. Code lists are also used to specify the address for data in tall-skinny format.

When mapping data using code lists to a form where multiple check boxes can be activated in one field, it is only possible to map the data based on choice number, not based on choice label. It is for this reason only possible to import values, not strings.

- Only the fields with a light blue background can be filled in.
- You can remove a table row by clicking on the trash can icon.
- You can import and edit an existing Define-XML file by clicking on **Import data mappings** in the Data mappings overview. Select the file you would like to import and click **Open**. Edit the data mapping table if necessary and click **Save changes** to save the data mapping.
- You can remove existing data mappings in the Data mappings overview by clicking **Delete**. It is not possible to delete a data mapping that has already been published.

#### Mapping the patient ID

5

- The subject ID should be mapped to {SubjectKey}, see image. Note that the mapping is case sensitive!
- Viedoc imports the data into existing subjects by matching the complete subject key as a string.
- If the data import should create new subjects, Viedoc needs to extract the country code, site code and subject sequence number from the subject ID. However, if the CSV file contains separate columns for country code or site code, Viedoc does not need to extract these from the subject ID.



• To extract the country code, site code and subject sequence number from the subject ID, it is necessary to map the format used for the subject ID, for example:

{CountryCode}-{SiteCode}-{SiteSubjectSeqNo}

The format of the subject ID is defined in the **Subject Id Generation Settings** in the **Study Settings** in Viedoc Designer. Click <u>here</u> for more information.

 Viedoc can only correctly extract the country code, site codes and subject sequence number within a subject ID if <u>one</u> of the following two requirements are met:

• The country code, site code and subject sequence number are separated with a separator, for example a hyphen: {CountryCode}-{SiteCode}-{SiteSubjectSeqNo}.

• The exact number of digits in the country code, the site code and the subject sequence number are specified in the subject ID generation settings in Viedoc Designer.

For example, if the subject ID generation settings in the study design are set at {CountryCode:00}{SiteCode:00} {SiteSubjectSeqNo:000}, the country code should consist of two digits, the site code should consist of two digits, and the subject sequence number should consist of three digits.

# 6 Mapping the visit ID and event date

- For scheduled events, the visit ID should be mapped to {StudyEventDefld}, see image. Note that the mapping is case sensitive!
- The date of the event can be mapped in two different ways:
  - Mapping to {EventDate}

{EventDate} is used to match the imported event date with the scheduled event, and to initiate events if the visit has not been initiated yet.



If the visit has been initiated, the Viedoc Import Application checks whether the date in the CSV file matches with the existing visit date. If the dates are matching, the data will be imported. If the dates are not matching, the data will not be imported.

#### • Mapping to {\$THIS.\$EVENT.EventDate}

{\$THIS.\$EVENT.EventDate} is used to update the event date. If the visit has been initiated, the Viedoc Import Application will update the existing date with the new date. Note that for matching the events, it is necessary to map the original event date to {EventDate} as well, so that the system recognizes which event should be updated.

- For unscheduled and common events, the event date should be mapped to {StudyEventDefld}, along with either {EventDate} or {StudyEventRepeatKey}. The event will then be matched either on event date or on event sequence number.
   If the event has not been initiated yet, the event date will be imported and/or an event sequence number will be created.
- You can also map data to a certain activity within a visit using the Activity ID. This is useful when the same form is used in two different activities within the same event, for example before and after administration of a drug. The Activity ID should be mapped to {FormRepeatKey}.

#### Mapping the data

- The data, for example laboratory results, should be mapped to the correct form and field ID, see image. The mapping of these data is also case sensitive.
- In the example in the image, the laboratory results of the alanine aminotransferase serum levels are mapped to {\$THIS.CC.RES\_ALAT}. This is build up as follows:

**\$THIS** maps the data to the visit as defined earlier on in the mapping.

CC is the ID of the form that the data should be imported into.

RES\_ALAT is the ID of the field in the form that the data should be imported into.

You can also specify a scheduled event explicitly in the data mapping. For example, it is possible to map the above data directly to Visit 1 using {V1.CC.RES\_ALAT}. In that case, the visit ID does not have to be included in the CSV file, and, hence, the visit ID does not need to be mapped to {StudyEventDefId}, because the event is already specified in the destination of the data.

D	A demo study									Save	changes	.0	
	domain   SAS domain name												
Cart	domain para domain name	,											
lata Napping Name						Domein neme							
900	tap example for Tall-Skiney3				e se	5 domair	rame						
	rted file structure						Viedoo						
	Column name	Description		Lick to			94	Destination (Medio) expression)			cr		
	PatientD	Subject ID	1										
		Subject ID	1					(Subject/kg)	1				
		Site Code Subject ID	1					KountryCodel-SiteCodel-SiteSubjectSepHol	1				
	visitio		1	None			۲	(SharlyEventDeffe)	1				
	VISITDAT	Voit Date	1				۲						
		Visit Date	1					[ExertDate]	1				
		Study Event Repeat Key	1					StudyEventRepeatSeyl	1				
	PARAMETER /		1	Nove			۲		1		۲		
	CLS ALAT							G1:	1	8			
	CLE AUP							0.2	1	1			
	CL3: AGAT							0.3	1	1			
	Retail 🖌		1	PARAMITER									
	ALAT	Alanine aminotransferase - Result	1					ISTHIS.CC.RES_ALATI	1				
2	ALB/P	Advaline phosphatase - Result						ISTHIS.CC.RES_ADAPI			۲		
3	·· 4547	Aspertate aminotransferase - Result	1					ISTHIS.CC.RES_ASAT)	1				
	Unit 🖌		1	PARAMETER									
1	ALAT	ALAT	1					(STHIS.CC.UNT_ALAT)	1				
2	·· 407	ADP	1					(STHS.CC.UNT_AUR)	1		۲		
3	ASAT	ASAT						ISTHIS.CC.UNT_AGATI					